# Premier Hb 9210™

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

Received June 2, 2014; accepted September 4, 2014

## Abstract

**Background:** The accurate and precise quantification of HbA<sub>1c</sub> is essential for the diagnosis and routine monitoring of patients with diabetes. We report an evaluation of the Trinity Biotech Premier Hb9210 analyser (Bray, Ireland/Kansas City, MO, USA), a boronate affinity chromatography-based high performance liquid chromatography (HPLC) system for the measurement of glycated haemoglobin.

which are well below the recommended upper limits of 3% CV for mmol/mol (SI units) and 2% CV for % (NGSP). The analyser showed a good correlation to HbA<sub>1c</sub> methods currently used in clinical practice and the IFCC reference method procedure. Haemoglobin variants AC, AS, AE and AD do not affect the measurement of HbA<sub>1c</sub>. Overall the Hb9210 performs well across the whole analytical range. **Conclusions:** The Hb9210 performs well and is suitable for clinical application in the analysis of HbA<sub>1c</sub>.

**Keywords:** affinity chromatography; diabetes; glycated haemoglobin; HbA<sub>1c</sub>; multicentre evaluation.

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# Manufacturing Plants

- **Acton, Massachusetts** – Supply raw materials to the research and diagnostic industries. 17 employees based at this site.
- **Buffalo, New York (IMMCO)** – Purchased in 2013, develops and manufactures autoimmune kits, complemented by specialized reference laboratory. 91 employees based at this site.
- **Jamestown, New York** – Specializes in production of Microtitre Plate EIA products for infectious diseases and auto-immunity. 116 employees based at this site.
- **Kansas City, Missouri** – Manufactures A1c range of products based on Boronate Affinity HPLC. 70 employees based at this site.
- **Carlsbad, California** – Specializes in development and manufacture of products utilizing Western Blot technology, including Lyme suite of products. 60 employees based at this site.
- **Belo Horizonte, Brazil** – Manufactures reagents for A1c range of products based on Boronate Affinity HPLC.
- **Uppsala, Sweden (Fiomi)** – Purchased in 2012, developing Troponin I and BNP tests for cardiac market. 26 employees based at this site.
- **Bray, Ireland** – Point-of-care / HIV, Immunofluorescence and Clinical Chemistry products manufactured at this site. 130 employees based at this site.

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**W. Garry John**

- Chairman of the IFCC Working Group
- Consultant and Honorary Professor of Clinical Biochemistry, Norwich & Norfolk University Hospital

**Randie Little**

- Coordinator of the NGSP Reference Network
- Research Professor, Co-Director of Diabetes Diagnostic Laboratory

**David B. Sacks**

- Chairman of the NGSP
- Chief of Clinical Chemistry, United States National Institute of Health

**Cas Weykamp**

- Coordinator of the IFCC Reference Network
- Director of MCA Laboratory, Queen Beatrix Hospital Winterswijk, NL

**Erna Lenters-Westra**

- Director of the IFCC Reference Laboratories
- HbA<sub>1c</sub> Research Coordinator, Isala, Zwolle, NL

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