


Automated Chemistry Analyzers



NormaChem 100 & 200

- Random access
- Benchtop, small footprint
- 130 resp. 200 tests per hour without ISE
- 200 resp. 275 tests per hour with ISE





Compact

Silent

Economic

NormaChem 100

NormaChem 200

TECHNICAL DATA

Throughput: 130 tests/hour w/o ISE
200 tests/hour with ISE

Sample volume: 2-35 μ L
28 sample positions

End-point, kinetic, immunoturbidimetry,
fixed-time, monochromatic, bi-chromatic,
multi standard analysis
Washing station
Non-stop cooling of reagent compartment
Reagent pre-heating

Throughput: 200 tests/hour w/o ISE
275 tests/hour with ISE

Sample volume: 3-30 μ L
44 sample positions

www.normadiagnostika.com

Norma Diagnostika GmbH

Ameisgasse 49-51/2A, 1140 Vienna, Austria • Phone: +43 1 914 8500 • Fax: +43 1 914 85 0715

Norma Diagnostika GmbH reserves the right without notice to make any technical changes it considers necessary or desirable.