



(°CO)

icon

Bernon

The Smart Hematology Analyzer

norme



- Microfluidic technology
- Patented, innovative solutions
- Trendsetting design

Will be launched soon.

Norma Icon has been designed for cost sensitive, environment conscious users

 (ϵ)

Smallest analyzer of its class
Smallest sample volume
Smallest reagent consumption

TECHNICAL DATA

Parameters: WBC, LYM, MID, GRA, LYM%, MID%, GRA%, HGB, RBC, HCT, MCV, RDWcv, RDWsd, MCH, MCHC, PLT, MPV, PCT, PDWcv, PDWsd, P-LCR%, P- LCC Histograms: WBC, RBC, PLT Sample volume: 9.6 μL whole blood (with EDTA) Sample tube: closed or open vials Throughput: 60 tests per hour Measuring technology: Volumetric impedance method and microfluidics Storage capacity: 100,000 test results incl. histograms User interface: 10.1 inch,1280 x 800 resolution LCD with touch screen Reagent consumption: 5.0 mL diluent, 1.0 mL lyse, 2.0 mL system solution

www.normadiagnostika.com

Norma Diagnostika GmbH

Hauptstr. 7., 3011 Untertullnerbach, Austria • Phone: +43 2233 52887 • Fax: +43 2233 52 88715

Norma Diagnostika GmbH reserves the right without notice to make any technical changes it considers necessary or desirable.