


norma
Diagnostic · Solutions

Electrolyte Analyzers




NormaLyte Series





NormaLyte Series

- Cost effective
- Low maintenance
- Low reagent consumption



TECHNICAL DATA

Technology: ISE (Ion Selective Electrode)

Sample: Serum / Plasma / Whole Blood / Urine / CSF

Sample Volume: 80-150 μ L

Reagent Volume: 1 mL

Throughput: 60 samples/hour

PARAMETERS

NormaLyte 100: K⁺, Na⁺

NormaLyte 200: K⁺, Na⁺, Cl

NormaLyte 300: K⁺, Na⁺, Cl, iCa

NormaLyte 400: K⁺, Na⁺, Cl, Li

NormaLyte 500: K⁺, Na⁺, Cl, iCa, TCa, pH

NormaLyte 600: K⁺, Na⁺, Cl, iCa, TCa, pH, TCO₂, Ag, KT

www.normadiagnostika.com

Norma Diagnostika GmbH

Ameisgasse 49-51/2A, 1140 Vienna, Austria • Phone: +43 1 914 8500 • Fax: +43 1 914 85 0715

Norma Diagnostika GmbH reserves the right without notice to make any technical changes it considers necessary or desirable.